

A 503

Mobile unit

No levels, for cleaning laboratory glassware/utensils.



- For loading modules and inserts
- Depending on the height, up to 6 modules can be used
- Upper rail and upper strut can be removed

EAN: 4002516854784 / Material number: 12685380 /  
Old Material Number: 69650302D

#### Product category

Large-capacity washer-disinfectors, laboratory

#### Machine affiliation

PLW 8615

PLW 8616

PLW 8617

#### Product category

Mobile unit with module connection

#### Application area

Reprocessing of laboratory glassware

#### Product characteristics

Material

Stainless steel, Plastic

Color

Stainless Steel, Black

#### Dimensions and weight

External dimensions, net height in inches

17 1/16 (433)

External dimensions, net width in inches

31 1/8 (790)

External dimensions, net depth in inches

25 5/8 (650)

External dimensions, gross height in inches

31 1/8 (790)

External dimensions, gross width in inches

28 3/16 (715)

External dimensions, gross depth in inches

33 7/8 (860)

Net weight in lbs

25 (11.40)

Gross weight in lbs

49 (22.20)