

Guard L1

canister vacuum cleaners

extremely powerful | ideal for carpets | intuitive operation



- ComfortFit: developed for the highest hygiene standards
- Durable: tested for the equivalent of 20 years** use
- Ideal power for every cleaning task: 4 selectable modes
- Individually selectable height and ergonomic Comfort handle
- Premium design concept with stylish applications
- Extra range, fewer socket changes: 36 ft (11 m) operating radius
- -
- AirClean filter: effective filtration for clean room air

* During the development phase of the advertised Guard vacuum cleaner series, Miele tested models and core components at various power levels for 750 hours of use (= 37.5 hours p.a.). More:

mieleusa.com/20years

EAN: 4002516855934 / Material number: 12704560 /

Old Material Number: 41UDE020CDN

General information	
Model number	Guard L1
Manufacturer material number	12704560
EAN number	4002516855934
Design	
Canister vacuum cleaner with dustbag	•
Product line	Guard L1
Motor type	PowerLine
Product type	SUDE0
Construction type and design	
Color	Nordic blue
User convenience	
Electronic suction power setting	+/- footswitches
Illuminated power indicator	•
Power indicator	Illuminated display
Telescopic suction wand	Suction wand X(HES)
Park system for vacuuming breaks	•
Park system for storing	on both sides
One-touch automatic cord rewind	•
Comfort handle 2.0	•
Mobility	
Caster wheels gentle on floors	•
4-point manoeuvrability	•
Operating radius in ft (m)	36 (11)
Cord length in ft (m)	24 (7.5)
Filter system	
FilterBag type	HyClean Pure TU
FilterBag™ volume in qt (US) (l)	3.7 (3.5) qt (US) (l)
FilterBag™ change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic FilterBag™ positioning	•
Technical data	
Power	1200
Volume in dB	74.0
Width in inches (mm)	19 5/8 (500)
Height in inches (mm)	9 3/8 (239)
Depth in inches (mm)	11 3/8 (288)
Weight in lbs (kg)	15 (6.5)
Standard accessories	
Upholstery nozzle	•
Crevice nozzle	•
Dust brush	•
Floorhead SBD 365-3	•
HyClean Air TU vacuum bag	•