

Guard L1 Cat & Dog canister vacuum cleaners powerful | attachment for removing pet hair | ideal for pet owners



- Active AirClean filter neutralizes unpleasant odors
- Durable: tested for the equivalent of 20 years** use
- ComfortFit: developed for the highest hygiene standards
- Individually selectable height and ergonomic Comfort handle
- For effortless cleaning: 2-in-1 upholstery and crevice nozzle
- Ideal power for every cleaning task: 4 selectable modes
- Premium design concept with stylish applications
- Turbobrush: removes stubborn soiling and animal hair

* During the development phase of the advertised Guard vacuum cleaner series, Miele tested models and core components at various power levels for 750 hours of use (= 37.5 hours p.a.). More: mieleusa.com/20years

EAN: 4002516855941 / Material number: 12704570 /
Old Material Number: 41UZE034CDN

Design	
Canister vacuum cleaner with dustbag	•
Product line	Guard L1
Motor type	PowerLine
Product type	SUZE0
Construction type and design	
Appliance color	Obsidian black
User convenience	
Electronic suction power setting	+/- footswitches
Illuminated power indicator	•
Telescopic suction wand	Suction wand X(HES)
Park system for vacuuming breaks	•
One-touch automatic cord rewind	•
Comfort handle 2.0	•
Mobility	
Caster wheels gentle on floors	•
4-point manoeuvrability	•
Operating radius in ft (m)	36 (11)
Cord length in ft (m)	24 (7.5)
Filter system	
FilterBag type	HyClean Pure TU
FilterBag™ volume in qt (US) (l)	3.7 (3.5) qt (US) (l)
FilterBag™ change indicator	•
Exhaust filter	Active AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic FilterBag™ positioning	•
Technical data	
Module depth in mm	28.80
Power	1200
Volume in dB	75.0
Width in inches (mm)	19 3/4 (500)
Height in inch (mm)	239
Weight in lbs (kg)	15 (6.6)
Standard accessories	
Accessory storage	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Dust brush	•
Floorhead SBD 365-3	•
Turbobrush STB 305-3 (DOUBLE; LOCKED)	•
HyClean Air TU vacuum bag	•