

APLW 908

Retrofitting kit

for dispensing cleaning agents from 10-liter media canisters.



- For using 10-liter media canisters outside the machine
- Suction lance for cleaning agents

EAN: 4002516930037 / Material number: 12755370 / Old Material Number: 69242020D

| Product category                                               |                |
|----------------------------------------------------------------|----------------|
| SlimLine laboratory glasswashers                               | •              |
| SlimLine laboratory glasswashers for cleaning industrial parts | •              |
| Machine affiliation                                            |                |
| PLW 7111                                                       | •              |
| PLW 7111 [Safety]                                              | •              |
| PWD 7121                                                       | •              |
| PWD 7122                                                       | •              |
| Product category                                               |                |
| Conversion kit                                                 | •              |
| Application area                                               |                |
| Reprocessing of laboratory glassware                           | •              |
| Reprocessing of dental instruments                             | •              |
| Product characteristics                                        |                |
| Material                                                       | Plastic        |
| Dimensions and weight                                          |                |
| External dimensions, gross height in inches                    | 3 3/4 (95)     |
| External dimensions, gross width in inches                     | 6 3/4 (170)    |
| External dimensions, gross depth in inches                     | 21 11/16 (550) |
| Net weight in lbs                                              | 1 (0.34)       |
| Gross weight in lbs                                            | 1 (0.50)       |