

ProCare Dent 10 A - 5 l
Liquid detergent, alkaline, 5 l
for the automatic reprocessing of dental instruments.



- Suitable for dental instruments, not for transmission instruments
- Adapted to the process sequences in Miele washer-disinfectors
- Excellent material/corrosion protection for stainless steel instruments
- Does not contain surfactants
- Free from oxidative additives

EAN: 4029577158194 / Material number: 12786490 / Old Material Number: 66912201V6

| | |
|---|-------------------------|
| Product category | |
| Thermal disinfectors | • |
| Applications | |
| Removal of blood | • |
| Removal of organic residues | • |
| Dispensing | |
| pH value at 68°F (20°C) | 14,0 |
| Dispensing recommendation | |
| Dispensing recommendation per rinse | 3 ml/l (0.3%) |
| Dispensing recommendation (per rinse) | 3 ml/l (0.3%) |
| Large opening for siphon | • |
| Intensity of soiling | Strong |
| Suitable for all levels of water hardness | • |
| Main applications | |
| Dental practices | • |
| Product characteristics | |
| Medium type | Alkaline cleaning agent |
| Aggregate form | Liquid concentrate |
| Ingredients | Phosphate, Silicates |
| Minimum storage temperature in °F/°C | -15 |
| Maximum storage temperature in °F/°C | 30 |
| Dimensions and weight | |
| External dimensions, net height in inches | 217 |
| External dimensions, net width in inches | 142 |
| External dimensions, net depth in inches | 248 |
| External dimensions, gross height in inches (mm) | 250.00 |
| External dimensions, gross width in inches (mm) | 190.00 |
| External dimensions, gross depth in inches (mm) | 150.00 |
| Gross weight in lbs (kg) | 9.64 |
| Fill volume in ml | 5000 |
| Product characteristics | |
| Medium type | Alkaline cleaning agent |
| Aggregate form | Liquid concentrate |
| Ingredients | Phosphate, Silicates |
| Minimum storage temperature in °F/°C | -15 |
| Maximum storage temperature in °F/°C | 30 |
| Suitable for automatic reprocessing of instruments and utensils | • |
| Suitable for standard instruments | • |
| Not suitable for transmission instruments | • |
| Does not contain oxidants | • |
| Does not contain surfactants | • |
| Good cleaning capability | • |
| High degree of material protection against corrosion | • |
| Dimensions and weight | |
| External dimensions, net height in inches | 8 9/16 (217) |
| External dimensions, net width in inches | 5 5/8 (142) |
| External dimensions, net depth in inches | 9 13/16 (248) |

ProCare Dent 10 A - 5 l
Liquid detergent, alkaline, 5 l
for the automatic reprocessing of dental instruments.



EAN: 4029577158194 / Material number: 12786490 / Old Material Number: 66912201V6

| | |
|---|--|
| External dimensions, gross height in inches | 9 7/8 (250) |
| External dimensions, gross width in inches | 7 1/2 (190) |
| External dimensions, gross depth in inches | 5 15/16 (150) |
| Gross weight in lbs | 0.35 (10) |
| Fill volume in ml | 169.07 (5000) fl oz (US) (Milliliter) |

| | |
|--|---|
| Compliance, test certificates, and symbols | |
| CE | • |
| GHS/CLP compliant | • |