

ProCare Dent 11 A - 10 kg  
 Powder cleaning agent, alkaline, 10 kg  
 for the automatic reprocessing of dental instruments.



- Suitable for dental instruments, not for transmission instruments
- Adapted to the process sequences in Miele washer-disinfectors
- Excellent material compatibility
- Free from oxidative additives
- Corrosion protection for stainless steel instruments and trays

EAN: 4029577158217 / Material number: 12786540 / Old Material Number: 66911201V6

<b>Product category</b>	
Thermal disinfectors	•
Large-capacity washer-disinfectors, dental	•
Washer-disinfectors, dental	•
<b>Applications</b>	
Removal of blood	•
Removal of organic residues	•
<b>Dispensing</b>	
Dispensing recommendation (per rinse)	3 g/l (30 g/cycle)
Manual dispensing	•
Intensity of soiling	Strong
Suitable for all levels of water hardness	•
<b>Main applications</b>	
Dental practices	•
<b>Product characteristics</b>	
Medium type	Alkaline cleaning agent
Aggregate form	Powder
Ingredients	Phosphate, Non-ionic surfactants, Silicates
pH value at 68°F (20°C)	12,0
Suitable for automatic reprocessing of instruments and utensils	•
Suitable for standard instruments	•
Not suitable for transmission instruments	•
Good cleaning capability	•
High degree of material protection against corrosion	•
<b>Storage information</b>	
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	30
<b>Dimensions and weight</b>	
External dimensions, net height in inches	10 9/16 (268)
External dimensions, net width in inches	10 7/16 (265)
External dimensions, net depth in inches	10 7/16 (265)
Load size in g	352.73 (10000)
Gross weight in lbs (kg)	23 (10)
<b>Compliance, test certificates, and symbols</b>	
CE	•
GHS/CLP compliant	•