

ProCare Dent 30 C - 5 l
Neutralizing agent, acidic, 5 l
for the automatic reprocessing of dental and transmission instruments.



- Suitable for dental instruments
- Adapted to the process sequences in Miele washer-disinfectors
- Good material compatibility, free from phosphates
- Neutralizes alkali residues that are carried over
- Based on citric acid

EAN: 4029577158224 / Material number: 12796330 /
Old Material Number: 66915303V6

Product category	
Thermal disinfectors	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing	
pH value at 68°F (20°C)	1,0
Dispensing recommendation	
Dispensing recommendation per rinse	1 ml/l (0.1%)
Dispensing recommendation (per rinse)	1 ml/l (0.1%)
Large opening for siphon	•
Suitable for all levels of water hardness	•
Main applications	
Dental practices	•
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Citric acid, Preservatives
Minimum storage temperature in °F/°C	-3
Maximum storage temperature in °F/°C	30
Dimensions and weight	
External dimensions, gross depth in inches (mm)	150.00
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Citric acid, Preservatives
Minimum storage temperature in °F/°C	-3
Maximum storage temperature in °F/°C	30
Suitable for anodized aluminum	•
Suitable for standard instruments and transmission instruments	•
Does not contain surfactants	•
Very good material compatibility	•
Suitable for sensitive materials	•
Neutralization of alkaline carryover	•
Excellent rinsing properties	•
Does not contain phosphorus compounds	•
Prevention of limescale formation in wash cabinet	•
Not suitable for transmission instruments	Suitable for standard/ transmission instruments, Suitable for anodized aluminum, Acid resistance Wash items and drains
Dimensions and weight	
External dimensions, net height in inches	8 9/16 (217)
External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)
External dimensions, gross height in inches	9 7/8 (250)
External dimensions, gross width in inches	7 1/2 (190)
External dimensions, gross depth in inches	5 15/16 (150)

ProCare Dent 30 C - 5 l
Neutralizing agent, acidic, 5 l
for the automatic reprocessing of dental and transmission instruments.



EAN: 4029577158224 / Material number: 12796330 /
Old Material Number: 66915303V6

Gross weight in lbs	0.21 (6)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)

Compliance, test certificates, and symbols	
CE	•
GHS/CLP compliant	•