

E 943/2

Module

for the optimum loading of narrow-neck glassware 100 - 500 ml.



- 32 nozzles (E 351) 4 x 160 mm with clips (E 353)

EAN: 4002516948452 / Material number: 12897120 / Old Material Number: 69594304D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert/module for narrow-neck glassware	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Capacity	
Narrow-neck glassware [number]	32