

E 944/2

Module

for the optimum loading of narrow-neck glassware 500 - 1,000 ml.



- 15 nozzles (E 352) 6 x 220 mm with clips (E 354)

EAN: 4002516948469 / Material number: 12897160 /
Old Material Number: 69594404D

Product category

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PG 8527

•

PG 8528

•

PLW 8615

•

PLW 8616

•

PLW 8617

•

Product category

Insert/module for narrow-neck glassware

•

Application area

Reprocessing of laboratory glassware

•

Product characteristics

Material

Stainless steel

Color

Stainless Steel

Capacity

Narrow-neck glassware [number]

15