

E 945/2 Module for inserts.



- Holding frame for inserts

EAN: 4002516949374 / Material number: 12897190 / Old Material Number: 69594504D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Other inserts/other modules	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel