

E 957

Mobile unit

for the optimum loading of up to 12 large volume laboratory glassware.

EAN: 4002516949695 / Material number: 12901310 / Old Material Number: 69595702D



- Equipped with 12 nozzles with Ø 10 mm
- Height adjustable holding frame
- Magnetic strip for automatic mobile unit recognition
- Max. vertical clearance from star overlay: 24 1/4 (615) in (mm) mm
- Requires A 857 conversion kit (not included) for
- use in PLW 861x

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Other mobile units	•
Module for large-volume laboratory glassware	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Maximum vertical clearance in mm	24 1/4 (615)
Capacity	
Large-volume laboratory glassware [number]	12