

ProCare Lab 10 AO - 5 l

liquid cleaning agent, alkaline, 5 L

for special machine cleaning of laboratory glassware and utensils.



- Note on the use of detergents with active chlorine *
- Oxidizing effect, free from surfactants
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics, and plastics
- Good cleaning for organic soiling

* When using detergents with active chlorine, increased wear of organic materials such as plastics and elastomers in the laboratory glassware washer is to be expected. During the main wash cycle with highly chlorinated detergents, a maximum temperature of 167°F (75°C) should not be exceeded. In order to avoid such elevated levels of wear and tear, Miele recommends that machines be serviced annually or after 1,000 operating hours.

A failure to perform regular maintenance may invalidate guarantee and warranty claims.

Please contact Miele Service or request an offer for a maintenance contract.

EAN: 4029577160265 / Material number: 12958890 /

Old Material Number: 62926102V2

Product category	
Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•
Applications	
Removal of protein	•
Removal of blood	•
Removal of residues of nutrient mediums	•
Removal of organic residues	•
Removal of cell and tissue culture residues	•
Removal of radioactive components	•
Removal of algae	•
Dispensing	
Recommended dosage per wash cycle	2-4 ml/l (0.2-0.4%)
Large opening for siphon	•
Automatic dispensing	•
Recommended water hardness level in mmol/L	0-1
Main applications	
Microbiological laboratories	•
Virological laboratories	•
Nuclear medicine laboratories	•
Laboratories for the food industry	•
Color industry	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Potassium hydroxide, Chlorine bleach, Phosphate
pH value at 68°F (20°C)	14,0
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Not suitable for aluminum	•
Not suitable for light and non-ferrous metals	•
Not suitable for light alloys	•
Free from microplastics	•
Oxidative effect due to active chlorine	•
Oxidizing effect	•
Storage information	
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	25
Dimensions and weight	
External dimensions, net height in inches	8 9/16 (217)
External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)
External dimensions, gross height in inches (mm)	9 7/8 (250)
External dimensions, gross width in inches (mm)	7 1/2 (190)
External dimensions, gross depth in inches (mm)	5 15/16 (150)

ProCare Lab 10 AO - 5 l
liquid cleaning agent, alkaline, 5 L
for special machine cleaning of laboratory glassware and utensils.



EAN: 4029577160265 / Material number: 12958890 / Old Material Number: 62926102V2	
Gross weight in lbs (kg)	15 (7.00)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
GHS/CLP compliant	•