

ProCare Lab 10 AP - 5 l

liquid cleaning agent, alkaline, 5 L

for optimum reprocessing of laboratory utensils.



- Alkaline cleaning agent with a high percentage of dispersants
- Free from surfactants, phosphates, and oxidants
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics, and plastics
- For cleaning stubborn or dried-on residues

EAN: 4029577160272 / Material number: 12958910 / Old Material Number: 62924102V2

Product category	
Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•
Applications	
Removal of grease	•
Removal of oil	•
Removal of organic residues	•
Dispensing	
Recommended dosage per wash cycle	2-4 ml/l (0.2-0.4%)
Large opening for siphon	•
Automatic dispensing	•
Recommended water hardness level in mmol/L	0-0.5
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral-oil industry	•
Industrial laboratories	•
Phosphate industry	•
Color industry	•
Laboratories for the pharmaceutical industry	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Potassium hydroxide, Polycarboxylates
pH value at 68°F (20°C)	14,0
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Not suitable for aluminum	•
Not suitable for light and non-ferrous metals	•
Does not contain oxidants	•
Does not contain phosphates	•
Does not contain surfactants	•
Free from microplastics	•
Good cleaning capability	•
High percentage of dispersants	•
Storage information	
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	30
Dimensions and weight	
External dimensions, net height in inches	8 9/16 (217)

ProCare Lab 10 AP - 5 l
liquid cleaning agent, alkaline, 5 L
for optimum reprocessing of laboratory utensils.



EAN: 4029577160272 / Material number: 12958910 / Old Material Number: 62924102V2

External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)
External dimensions, gross height in inches (mm)	9 7/8 (250)
External dimensions, gross width in inches (mm)	7 1/2 (190)
External dimensions, gross depth in inches (mm)	5 15/16 (150)
Gross weight in lbs (kg)	17 (7.50)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)

Compliance, test certificates, and symbols	
GHS/CLP compliant	•