

ProCare Lab 10 AP - 10 l

liquid cleaning agent, alkaline, 10 L

for universal machine cleaning of laboratory glassware and utensils.



- Alkaline cleaning agent with a high percentage of dispersants
- Free from surfactants, phosphates, and oxidants
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics, and plastics
- For cleaning stubborn or dried-on residues

EAN: 4029577160289 / Material number: 12958930 /

Old Material Number: 62924202V2

Product category	
Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•
Applications	
Removal of grease	•
Removal of oil	•
Removal of organic residues	•
Dispensing	
Recommended dosage per wash cycle	2-4 ml/l (0.2-0.4%)
Large opening for siphon	•
Automatic dispensing	•
Recommended water hardness level in mmol/L	0-0.5
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral-oil industry	•
Industrial laboratories	•
Phosphate industry	•
Color industry	•
Laboratories for the pharmaceutical industry	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Potassium hydroxide, Polycarboxylates
pH value at 68°F (20°C)	14,0
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Not suitable for aluminum	•
Not suitable for light and non-ferrous metals	•
Not suitable for light alloys	•
Does not contain oxidants	•
Does not contain phosphates	•
Does not contain surfactants	•
Free from microplastics	•
Good cleaning capability	•
High percentage of dispersants	•
Storage information	
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	30

ProCare Lab 10 AP - 10 l
liquid cleaning agent, alkaline, 10 L
for universal machine cleaning of laboratory glassware and utensils.



EAN: 4029577160289 / Material number: 12958930 /
Old Material Number: 62924202V2

Dimensions and weight	
External dimensions, net height in inches	12 3/16 (308)
External dimensions, net width in inches	7 9/16 (191)
External dimensions, net depth in inches	9 1/4 (234)
External dimensions, gross height in inches (mm)	12 5/8 (320)
External dimensions, gross width in inches (mm)	9 1/16 (230)
External dimensions, gross depth in inches (mm)	7 1/2 (190)
Gross weight in lbs (kg)	30 (13.50)
Fill volume in ml	338.15 (10000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
GHS/CLP compliant	•