

ProCare Lab 11 AP - 10 kg

Powder cleaning agent, alkaline, 10 kg

for universal machine cleaning of laboratory glassware and utensils.



- Alkaline cleaning agent with oxidizing effect
- Free from phosphates, contains oxygen-based bleaching agents
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics, and plastics
- Powder cleaning agent for manual dispensing

EAN: 4029577160326 / Material number: 12960520 /

Old Material Number: 62922102V2

Product category	
Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•
Applications	
Removal of blood	•
Removal of organic residues	•
Dispensing	
Recommended dosage per wash cycle	2-4 g/l (20-40 g/cycle)
Manual dispensing	•
Recommended water hardness level in mmol/L	0-1.1
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral-oil industry	•
Serological institutes	•
Phosphate industry	•
Color industry	•
Laboratories for the pharmaceutical industry	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Powder
Ingredients	Oxygen-based bleaching agent, Surfactants, Silicates
pH value at 68°F (20°C)	13,0-14,0
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Not suitable for aluminum	•
Not suitable for light and non-ferrous metals	•
Not suitable for light alloys	•
Does not contain phosphates	•
Free from microplastics	•
Good cleaning capability	•
Material protection	•
Oxidizing effect	•
Storage information	
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	25
Dimensions and weight	
External dimensions, net height in inches	10 9/16 (268)
External dimensions, net width in inches	10 7/16 (265)
External dimensions, net depth in inches	10 7/16 (265)

ProCare Lab 11 AP - 10 kg

Powder cleaning agent, alkaline, 10 kg

for universal machine cleaning of laboratory glassware and utensils.



EAN: 4029577160326 / Material number: 12960520 /

Old Material Number: 62922102V2

External dimensions, gross height in inches (mm)	11 1/16 (280)
External dimensions, gross width in inches (mm)	10 11/16 (270)
External dimensions, gross depth in inches (mm)	10 11/16 (270)
Load size in g	352.73 (10000)
Gross weight in lbs (kg)	23 (10.50)

Compliance, test certificates, and symbols

GHS/CLP compliant

•