

ProCare Lab 30 C - 5 l

Neutralizing agent based on citric acid, acidic, 5 L

for material-saving machine cleaning of laboratory glassware and utensils.



- To be used as neutralizing agent or acid pre-cleaning agent
- Free from surfactants and phosphates
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics, and plastics
- Removes residues that are soluble in acids from sensitive materials

EAN: 4029577160340 / Material number: 12960640 /
Old Material Number: 62927102V2

Product category

Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•

Applications

Removal of limescale and other acid-soluble impurities	•
--	---

Dispensing

Recommended dosage per wash cycle	1–5 ml/l (0.1–0.5 %)
Large opening for siphon	•
Automatic dispensing	•
Suitable for all levels of water hardness	•

Main applications

Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the phosphate industry	•
Laboratories for the pharmaceutical industry	•

Product characteristics

Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Citric acid, Preservatives
pH value at 68°F (20°C)	1,0
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Suitable for sensitive materials	•
Suitable for neutralization of non-anodized aluminum	•
Not suitable for anodized aluminum	•
Not suitable for light alloys	•
Not suitable for acidic pre-cleaning of non-anodized aluminum	•
Does not contain phosphates	•
Does not contain surfactants	•
Free from microplastics	•
Very gentle on materials	•
Acidic pre-cleaning	•
Neutralization of alkaline carryover	•

Storage information

Minimum storage temperature in °F/°C	-3
Maximum storage temperature in °F/°C	30

Dimensions and weight

External dimensions, net height in inches	8 9/16 (217)
External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)

ProCare Lab 30 C - 5 l
Neutralizing agent based on citric acid, acidic, 5 L
for material-saving machine cleaning of laboratory glassware and utensils.



EAN: 4029577160340 / Material number: 12960640 / Old Material Number: 62927102V2	
External dimensions, gross height in inches (mm)	9 7/8 (250)
External dimensions, gross width in inches (mm)	7 1/2 (190)
External dimensions, gross depth in inches (mm)	5 15/16 (150)
Gross weight in lbs (kg)	14 (6.50)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
GHS/CLP compliant	•