

Guard L1 Allergy canister vacuum cleaners powerful | HEPA AirClean filter | ideal for allergy sufferers



- -
- Ideal power for every cleaning task: 4 selectable modes
- Premium design concept with stylish applications
- Individually selectable height and ergonomic Comfort handle
- For effortless cleaning: 2-in-1 upholstery and crevice nozzle
- Ideal for allergy sufferers thanks to HEPA AirClean filter
- -
- ComfortFit: developed for the highest hygiene standards

EAN: 4002517026203 / Material number: 13115840 /
Old Material Number: 41UZE033CDN

General information	
Model number	Guard L1 Allergy
Manufacturer material number	13115840
EAN number	4002517026203
Design	
Canister vacuum cleaner with dustbag	•
Product line	Guard L1
Product type	SUZE0
Construction type and design	
Color	Brilliant White
User convenience	
Electronic suction power setting	+/- footswitches
Power indicator	Illuminated display
Telescopic suction wand	Suction wand X(HES)
Park system for vacuuming breaks	•
Park system for storing	on both sides
One-touch automatic cord rewind	•
Mobility	
Caster wheels gentle on floors	•
Three-point maneuverability	•
Operating radius in ft (m)	39 (12)
Cord length in ft (m)	27 (8.5)
Filter system	
FilterBag type	HyClean Pure TU
FilterBag™ volume in qt (US) (l)	3.7 (3.5) qt (US) (l)
FilterBag™ change indicator	•
Exhaust filter	HEPA AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic FilterBag™ positioning	•
Technical data	
Power	890
Volume in dB	75.0
Width in inches (mm)	19 5/8 (500)
Height in inches (mm)	9 3/8 (239)
Depth in inches (mm)	11 3/8 (288)
Weight in lbs (kg)	15 (6.6)
Standard accessories	
Upholstery nozzle	•
Crevice nozzle	•
Dust brush	•
Floorhead SBD 365-3	•
HyClean Air TU vacuum bag	•
Voucher for HyClean Air TU	•