


E 328

Support rail

for the optimum loading of upright instruments.



- The insert can hold upright instruments

EAN: 4002511117051 / Material number: 02795940 / Old Material Number: 69732801D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
<b>Product category</b>	
Other accessories	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	15
External dimensions, net width in mm	34
External dimensions, net depth in mm	145
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0,06
Gross weight in kg	0,07