

E 379

Insert

For holding different instruments.



- 0.8 mm wire mesh, 1.7 mm mesh width, 5 mm frames
- equipped with 2 handles
- For use in upper and lower basket

EAN: 4002511116450 / Material number: 03203500 / Old Material Number: 69537901D

| Product category | |
|---|-----------------|
| Thermal disinfectors | • |
| Washer-disinfectors, medical | • |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7826 | • |
| PG 8527 | • |
| PG 8528 | • |
| G 7831 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8581 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8591 | • |
| PG 8592 | • |
| Product category | |
| Mesh tray | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of ENT instruments | • |
| Reprocessing of dental instruments | • |
| Reprocessing of surgical instruments | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0,00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 110 |
| External dimensions, net width in mm | 180 |
| External dimensions, net depth in mm | 445 |
| External dimensions, gross height in mm | 135 |
| External dimensions, gross width in mm | 235 |
| External dimensions, gross depth in mm | 485 |
| Net weight in kg | 0,72 |
| Gross weight in kg | 1,25 |