

E 146

Insert

For the optimum loading of different instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 8 mm
- equipped with 2 hinged handles
- For use in upper and lower basket

EAN: 4002511116009 / Material number: 03346680 / Old Material Number: 69514601D

<b>Product classification</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Mesh tray	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
Reprocessing of surgical instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	55
External dimensions, net width in mm	150
External dimensions, net depth in mm	225
External dimensions, gross height in mm	75
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	245
Net weight in kg	0,61
Gross weight in kg	0,70