

E 135/2

Insert

for optimum loading of 19 baby bottles up to 90 ml.



- Suitable for 19 baby bottles x 90 ml
- Suitable for bottles sized 51 x 51 mm, neck: 42 x 42 mm
- For upper basket A 102 and lower basket A 151

EAN: 4002511115927 / Material number: 03808650 / Old Material Number: 69513503D

| | |
|---|-----------------|
| Product category | |
| Washer-disinfectors, medical | • |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7823 | • |
| G 7824 | • |
| G 7826 | • |
| PG 8528 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8592 | • |
| Product category | |
| Insert for various flasks | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of baby bottles | • |
| Capacity | |
| Baby bottles [number] | 19 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0,00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 125 |
| External dimensions, net width in mm | 192 |
| External dimensions, net depth in mm | 447 |
| External dimensions, gross height in mm | 140 |
| External dimensions, gross width in mm | 235 |
| External dimensions, gross depth in mm | 480 |
| Net weight in kg | 1,97 |
| Gross weight in kg | 2,10 |