

E 329

Mobile unit

for the optimum loading of narrow necked glassware.



- equipped with 39 injector nozzles

EAN: 400251116108 / Material number: 03809340 / Old Material Number: 69532901D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Mobile unit for narrow-necked glassware	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	39
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	517
External dimensions, gross height in mm	545
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	530
Net weight in kg	5,50
Gross weight in kg	6,20
<b>Companion pack</b>	
ID 90, injector nozzle, Ø 2.5 mm, length 90 mm [number]	4
ID 110, injector nozzle, Ø 2.5 mm, length 110 mm [number]	5
ID 140, injector nozzle, Ø 4 mm, length 140 mm [number]	5
ID 160, injector nozzle, Ø 4 mm, length 160 mm [number]	5
ID 180, injector nozzle, Ø 4 mm, length 180 mm [number]	5
ID 200, injector nozzle, Ø 6 mm, length 200 mm [number]	5
ID 220, injector nozzle, Ø 6 mm, length 220 mm [number]	5
ID 220, injector nozzle, Ø 6 mm, length 240 mm [number]	5
Rinse nozzle for powder dispenser [number]	1