

E 336

Irrigation sleeves

For the optimum loading of pipettes/MIS instruments in injector units.

EAN: 4002511117068 / Material number: 03809390 / Old Material Number: 69733601D



- Pipettes up to max. 445 mm length can be held
- Plastic, screw connection

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Irrigation sleeve	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
<b>Capacity</b>	
Pipettes [number]	1
Lumened instruments [number]	1
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	121
External dimensions, net width in mm	11
External dimensions, net depth in mm	11
External dimensions, gross height in mm	10
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	125
Diameter in mm	11
Net weight in kg	0,02
Gross weight in kg	0,02