

E 367

Mobile unit

for the optimum loading of intubation material.



- 20 nozzles Ø 2.5 x 30 mm, 25 nozzles Ø 4.0 x 30 mm (5 loose)
- 40 holding springs for nozzles
- Equipped with connection for hot air drying
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002511116337 / Material number: 03809710 / Old Material Number: 69536701D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit for anaesthetic instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of anaesthetic instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	515
External dimensions, gross height in mm	490
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	550
Net weight in kg	6,20
Gross weight in kg	10,00