

E 406
Mobile unit
for the optimum loading of up to 116 pipettes.



- Suitable for 116 pipettes up to 450 mm
- Surface area 16 x 16 mm
- Complete with frame height 150 mm
- Without bracket for magnetic strip for mobile unit recognition

EAN: 4002511116665 / Material number: 03810110 / Old Material Number: 69540601D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit for pipettes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Pipettes [number]	116
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	517
External dimensions, gross height in mm	500
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	530
Net weight in kg	4,90
Gross weight in kg	6,90