


ID 240  
 Injector nozzle  
 with plastic supports, Ø 6 mm, length 240 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring flasks and test tubes, for example

EAN: 4002511117037 / Material number: 03810400 / Old Material Number: 69724001D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	55
External dimensions, net width in mm	55
External dimensions, net depth in mm	240
External dimensions, gross height in mm	60
External dimensions, gross width in mm	60
External dimensions, gross depth in mm	240
Diameter in mm	6
Net weight in kg	0,03
Gross weight in kg	0,04