

E 416

Insert

for optimum loading of 6 single and two-piece specula.



- Suitable for 6 single and two-part specula
- 7 holders 40 mm apart
- For use in upper (only with single-part specula) and lower basket

EAN: 4002511803824 / Material number: 03888650 / Old Material Number: 69541601D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
<b>Product category</b>	
Insert for gynaecological instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of GYN instruments	•
<b>Capacity</b>	
GYN specula [number]	6
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	157
External dimensions, net width in mm	178
External dimensions, net depth in mm	279
External dimensions, gross height in mm	210
External dimensions, gross width in mm	200
External dimensions, gross depth in mm	365
Net weight in kg	0,69
Gross weight in kg	0,97