

E 442

Irrigation sleeves

for the optimum cleaning of MIS instruments.



- Suitable for instruments with Ø 4 – 8 mm, screw connection
- Length 121 mm mm, Ø 6 mm mm
- Cap, opening Ø 6 mm
- Mounting clip (mat. no. 4174850)

EAN: 4002511908208 / Material number: 04203830 / Old Material Number: 69744201D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Irrigation sleeve	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	121
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Diameter in mm	6
Net weight in kg	0,03
Gross weight in kg	0,03