

E 446

Irrigation sleeve cover

for the optimum securing of instruments in irrigation sleeves.



- 12 caps for irrigation sleeves
- Opening Ø 10 mm mm

EAN: 4002511908178 / Material number: 04203860 / Old Material Number: 69744601D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Other accessories	•
Product properties	
Material	Plastic
Colour	Transparent
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Diameter in mm	10
Net weight in kg	0,03
Gross weight in kg	0,03