

E 464

Holder

For E 454 injector nozzles Ø 13 x 65 mm.



- Equipped with a spring for height adjustment

EAN: 4002512173018 / Material number: 04837600 / Old Material Number: 69746401D

| Product classification | |
|---|-----------------|
| Washer-disinfectors, medical | • |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7823 | • |
| G 7824 | • |
| G 7825 | • |
| G 7826 | • |
| PG 8527 | • |
| PG 8528 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8592 | • |
| Product category | |
| Holder | • |
| Injector nozzle | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of MIS instruments | • |
| Reprocessing of lumened instruments | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0,00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 65 |
| External dimensions, net width in mm | 0 |
| External dimensions, net depth in mm | 0 |
| External dimensions, gross height in mm | 12 |
| External dimensions, gross width in mm | 80 |
| External dimensions, gross depth in mm | 120 |
| Diameter in mm | 13 |
| Net weight in kg | 0,04 |
| Gross weight in kg | 0,05 |