

E 930

Insert

for optimum loading of up to 36 theatre shoes up to size 44.



- For use in lower level of mobile unit E 975/2
- Equipped with 36 holders, height 280 mm

EAN: 4002512257145 / Material number: 05001360 / Old Material Number: 69593001D

<b>Product classification</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Insert for theatre shoes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Capacity</b>	
Theatre shoes [number]	36
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	289
External dimensions, net width in mm	586
External dimensions, net depth in mm	752
External dimensions, gross height in mm	350
External dimensions, gross width in mm	770
External dimensions, gross depth in mm	1000
Net weight in kg	3,65
Gross weight in kg	5,70