

## AUF 1 Holder

for optimum reprocessing of transmission instruments in upper basket O 177/1.



- Comprises: seal, insert nut, clamp
- Rack lower part, rack upper part, 5 filter plates

EAN: 4002512427517 / Material number: 05159510 / Old Material Number: 68551101D

<b>Product category</b>	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8591	•
<b>Product category</b>	
Holder and/or adapter for transmission instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	80
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	100
Net weight in kg	0,90
Gross weight in kg	1,84
<b>Companion pack</b>	
FP filter plate	•