

E 713

Mobile unit

for loading of 2 DIN containers including lids as well as 4 mesh trays on 3 levels.



- Equipped with 2 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (without magnets ML)
- Equipped with a connection for 5 hollow instruments
- Equipped with connection for hot air drying

EAN: 4002512401944 / Material number: 05221550 / Old Material Number: 69571301D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7825	•
G 7826	•
Product category	
Mobile unit for inserts	•
Mobile unit for containers	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of surgical instruments	•
Reprocessing of containers	•
Capacity	
DIN mesh trays [number]	4
SPRI mesh trays [number]	2
DIN containers, height 15 cm [number]	2
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	530
External dimensions, net width in mm	575
External dimensions, net depth in mm	600
External dimensions, gross height in mm	560
External dimensions, gross width in mm	640
External dimensions, gross depth in mm	680
Net weight in kg	9,70
Gross weight in kg	10,00
Companion pack	
Magnetic strip for automatic mobile unit recognition (without magnets) [number]	1