


E 475  
Holding bar  
for stabilising sleeves.



- Suitable for MIS inserts

EAN: 4002512401234 / Material number: 05226040 / Old Material Number: 69747501D

<b>Product classification</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Holder	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	110
External dimensions, net width in mm	155
External dimensions, net depth in mm	500
External dimensions, gross height in mm	110
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	510
Net weight in kg	0,10
Gross weight in kg	0,22