


E 489

Universal bracket

for loading diverse wash items such as insoles for theatre shoes.



EAN: 4002512431835 / Material number: 05337120 / Old Material Number: 69548901D

| Product classification | |
|---|-----------------|
| Washer-disinfectors, medical | • |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7823 | • |
| G 7824 | • |
| G 7825 | • |
| G 7826 | • |
| PG 8528 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8592 | • |
| Product category | |
| Insert for theatre shoes | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of theatre shoes | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0,00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 60 |
| External dimensions, net width in mm | 464 |
| External dimensions, net depth in mm | 10 |
| External dimensions, gross height in mm | 35 |
| External dimensions, gross width in mm | 130 |
| External dimensions, gross depth in mm | 470 |
| Net weight in kg | 0,12 |
| Gross weight in kg | 0,23 |