


E 984  
Insert  
for a wide variety of utensils.



- Can be equipped with different holders

EAN: 4002512431323 / Material number: 05337750 / Old Material Number: 69598401D

<b>Product classification</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Insert for kidney dishes/basins	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	371
External dimensions, net width in mm	150
External dimensions, net depth in mm	572
External dimensions, gross height in mm	390
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	590
Net weight in kg	2,95
Gross weight in kg	3,30