


E 989

Universal bracket

for 8 insoles for theatre shoes, for example.



- For use in insert E 984

EAN: 4002512431828 / Material number: 05376740 / Old Material Number: 69598901D

Product classification	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8528	•
Product category	
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of theatre shoes	•
Capacity	
Insoles [number]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	20
External dimensions, gross width in mm	330
External dimensions, gross depth in mm	415
Net weight in kg	0,10
Gross weight in kg	0,26