


A 18  
Lid  
For fixing delicate wash items.



- 1 mm wire, 5 mm mesh size, 4 mm frames

EAN: 4002513118124 / Material number: 05631080 / Old Material Number: 69745006D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Lid	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	20
External dimensions, net width in mm	230
External dimensions, net depth in mm	430
External dimensions, gross height in mm	70
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	490
Net weight in kg	0,37
Gross weight in kg	0,55