

E 477
 Holders
 for micro instruments.



- For use in trays with 5 mm mesh
- 20 items per bag

EAN: 4002513239119 / Material number: 05692390 / Old Material Number: 69747701D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other accessories	•
Product properties	
Material	Plastic
Colour	Grey
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0,03
Gross weight in kg	0,04