

E 492
Insert
for the optimum loading of kidney trays.



- equipped with 9 holders, height 86 mm mm, distance 49 mm mm
- For use in lower basket

EAN: 4002513317343 / Material number: 05771140 / Old Material Number: 69549201D

| | |
|---|-----------------|
| Product category | |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7831 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8592 | • |
| Product category | |
| Insert for kidney dishes/basins | • |
| Other inserts | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Distance between rails in mm | 49 |
| Rail height in mm | 86 |
| Area of application | |
| Reprocessing of kidney dishes | • |
| Capacity | |
| Kidney dishes [number] | 9 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0,00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 120 |
| External dimensions, net width in mm | 256 |
| External dimensions, net depth in mm | 474 |
| External dimensions, gross height in mm | 125 |
| External dimensions, gross width in mm | 320 |
| External dimensions, gross depth in mm | 480 |
| Net weight in kg | 1,06 |
| Gross weight in kg | 1,31 |