

E 710/1

Mobile unit

for optimum loading of 4 DIN containers including lids.



- Container size: 150 x 300 x 600 mm
- Lid size: 30 x 300 x 600 mm
- Magnetic strip for mobile unit recognition (without magnets ML)

EAN: 4002513369199 / Material number: 05848260 / Old Material Number: 69571002D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Mobile unit for containers	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of containers	•
<b>Capacity</b>	
DIN containers, height 15 cm [number]	4
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	450
External dimensions, net width in mm	530
External dimensions, net depth in mm	600
External dimensions, gross height in mm	740
External dimensions, gross width in mm	800
External dimensions, gross depth in mm	1200
Net weight in kg	29,00
Gross weight in kg	38,00
<b>Companion pack</b>	
Magnetic strip for automatic mobile unit recognition (without magnets) [number]	1