

UE 30-30/60-78

Box plinth

to raise the side unit, 60 cm deep, 30 cm wide.



- Suitable for side cabinets G 7895 and G 7896

EAN: 4002513369861 / Material number: 05848340 / Old Material Number: 69373004D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Product category	
Plinth/base	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	300
External dimensions, net width in mm	300
External dimensions, net depth in mm	600
External dimensions, gross height in mm	380
External dimensions, gross width in mm	400
External dimensions, gross depth in mm	630
Net weight in kg	7,00
Gross weight in kg	8,10