

E 441/1

Insert

for the optimum loading of Micro instruments.



- The mesh width on the base is 1.7 mm
- The sides are closed, stackable
- Sections w/ 6 adjustable dividers for positioning of instruments
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002513460476 / Material number: 05988240 / Old Material Number: 69544102D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert/module for MIS instruments	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of MIS instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	60
External dimensions, net width in mm	183
External dimensions, net depth in mm	284
External dimensions, gross height in mm	100
External dimensions, gross width in mm	210
External dimensions, gross depth in mm	315
Net weight in kg	0,77
Gross weight in kg	0,92