

MF/3
Trolleys
for loading & unloading washer-&disinfectors.



- For easy handling of the mobile instrument unit
- With lifting technology operated by foot pressure
- Equipped with 2 lockable castors
- Infeed height 751 mm, +/- 100 mm

EAN: 4002513771114 / Material number: 06392900 / Old Material Number: 69200107D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
Product category	
Storage/transport system	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	1182
External dimensions, net width in mm	660
External dimensions, net depth in mm	807
External dimensions, gross height in mm	430
External dimensions, gross width in mm	800
External dimensions, gross depth in mm	1200
Net weight in kg	34,20
Gross weight in kg	52,40