


U 874/1

Lower basket

for the optimum loading of different inserts and mesh trays.



- Vertical clearance can vary with different upper baskets
- With bracket for magnetic strip for ML/2 autom. unit recognition

EAN: 4002513808919 / Material number: 06442940 / Old Material Number: 69187403D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Carrier/basket carrier	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of MIS instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	50
External dimensions, net width in mm	531
External dimensions, net depth in mm	515
External dimensions, gross height in mm	100
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	545
Net weight in kg	2,10
Gross weight in kg	2,97