

E 729

Insert

for optimum loading of up to 4-5 ophthalmic theatre sets on 3 levels.



- 2 levels for 4 DIN mesh trays
- 1 level with 16 or 30 injector nozzles
- Connection option for E 478/1 to the inlet pipe
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (without magnets ML)

EAN: 4002513893274 / Material number: 06517580 / Old Material Number: 69572901D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Mobile unit for ophthalmology instruments	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of ophthalmology instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	4
SPRI mesh trays [number]	2
Theatre sets for micro instruments [number]	4-5
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	471
External dimensions, net width in mm	530
External dimensions, net depth in mm	600
External dimensions, gross height in mm	680
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	650
Net weight in kg	14,40
Gross weight in kg	16,40
<b>Companion pack</b>	
E 476 holder [number]	1
E 477 stopper [number]	1
E 790 Luer Lock connection element female to female [number]	8
E 791 Luer Lock connection element female to male [number]	8
Magnetic strip for automatic mobile unit recognition (without magnets) [number]	1