

E 901/2

Mobile unit

for the optimum loading of mesh trays/basins and hollow instruments.



- Equipped with 3 integrated spray arms
- The loading areas for the individual levels are different
- The loading capacity for the individual levels are different
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (incl. magnets ML)

EAN: 4002514219882 / Material number: 06893870 / Old Material Number: 69590103D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for MIS instruments	•
Mobile unit for basins	•
Mobile unit for mesh trays	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Capacity	
DIN mesh trays [number]	12
SPRI mesh trays [number]	4
Lumened instruments [number]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	590
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	720
External dimensions, gross width in mm	830
External dimensions, gross depth in mm	780
Net weight in kg	26,50
Gross weight in kg	30,00
Companion pack	
E 362 blind stopper [number]	3
E 447 female adapter, for male Luer Lock [number]	2
E 449 male adapter, for female Luer Lock [number]	2
E 980 angled adapter [number]	3
Injector nozzle, Ø 4 mm, length 110 mm [number]	4
Injector nozzle, Ø 4 mm, length 120 mm, with spring clip [number]	3
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1