

E 913/3

Mobile unit

for the optimum loading of containers, lids and mesh trays on 3 levels.



- The loading areas for the individual levels are different
- Equipped with 2 integrated spray arms
- Equipped with connections for 5 hollow instruments
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (incl. magnets ML)

EAN: 4002514219899 / Material number: 06893880 / Old Material Number: 69591304D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for containers	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of surgical instruments	•
Reprocessing of containers	•
Capacity	
DIN mesh trays [number]	6
SPRI mesh trays [number]	2
DIN containers, height 15 cm [number]	3
Lumened instruments [number]	5
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	574
External dimensions, net width in mm	790
External dimensions, net depth in mm	640
External dimensions, gross height in mm	720
External dimensions, gross width in mm	830
External dimensions, gross depth in mm	780
Net weight in kg	26,50
Gross weight in kg	30,00
Companion pack	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1