


E 906/1

Module

for optimum loading of MIS instruments, such as arthroscopy and urology.



- Holds up to 10 hollow instruments
- Subdivision for arthroscopy and laparoscopy, for example.

EAN: 4002514256313 / Material number: 06943250 / Old Material Number: 69590602D

| | |
|--|-----------------|
| Product category | |
| Washer-disinfectors, medical | • |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7823 | • |
| G 7824 | • |
| G 7825 | • |
| G 7826 | • |
| PG 8527 | • |
| PG 8528 | • |
| PG 8536 | • |
| Product category | |
| Insert/module for MIS instruments | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of MIS instruments | • |
| Capacity | |
| Lumened instruments [number] | 10 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0,00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 40 |
| External dimensions, net width in mm | 461 |
| External dimensions, net depth in mm | 510 |
| External dimensions, gross height in mm | 160 |
| External dimensions, gross width in mm | 545 |
| External dimensions, gross depth in mm | 530 |
| Net weight in kg | 2,87 |
| Gross weight in kg | 3,56 |
| Companion pack | |
| E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number] | 1 |
| E 362 blind stopper [number] | 2 |
| E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number] | 5 |
| E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number] | 3 |
| E 448 silicone hose with female Luer Lock adapter [number] | 2 |
| E 454 injector nozzle for trocar sleeves 10–15 mm [number] | 1 |
| E 456 spring stay for MIS instruments [number] | 2 |
| E 464 holder for injector nozzle E 454 [number] | 1 |
| E 908/1 insert for separable MIS instruments/work inserts [number] | 1 |