

E 949

Module

for the optimum loading of 8 intramedullary nails and/or medullary reamers.



- Equipped with 8 holders for intramedullary nails*
- Equipped with 6 extension sleeves for medullary reamers**

* The holders are pre-installed.

** The extension sleeves are supplied loose.

EAN: 4002514312040 / Material number: 07025320 / Old Material Number: 69594901D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of orthopaedic surgery instruments	•
Capacity	
Bone marrow drills [number]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	115
External dimensions, net width in mm	352
External dimensions, net depth in mm	558
External dimensions, gross height in mm	290
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	560
Net weight in kg	35,00
Gross weight in kg	45,00